

Areas of Disinfection	Sub-Area	Device	Level of Disinfection	Frequency	Who	Disinfection Agent
<b>Professional</b>	Exam Room	Gonioscopy Lenses, Tonometer/Pachimeter probes	High	3	A	2
		Contact Lenses	Standard	1	A	3
		Forehead/Chin rests (phoropter, perimeter, OCT, camera, auto-refractor)	High	1	A, B	1
		Occuluders, Eye patches	Standard	1	A	1
		Diagnostic Equipment (i.e perimeter , OCT)	Standard	1	A, B	1
		Sinks	Standard	2	B	4
		Exam Chair and Unit	Standard	2	AB	4
		Trial Frame	Standard	3	A	1
		Handheld Instruments	Standard	3	A	1
		Lab/Dispensing Area	Contact Lens Cases	High	1	A, B
	Frame warmer		Standard	4	B	1
	Frames on Display		Standard	2	B	1
	Lab Hand Tools		Standard	4	B	1
	<b>Administrative</b>		Desk Counters	Standard	2	B
Computer Keyboards, Mouse & Telephone			Standard	4	B	1, 4
VISA Device			Standard	2	B	1, 4
Staplers, Tape Dispensers			Standard	4	B	1, 4
Pens, Pencils			Standard	2	B	1
Fax Machine			Standard	4	B	1, 4
Waiting Area			Standard	2	B	1, 4
Door Handles			Standard	2	B	1, 4
Light Switches			Standard	2	B	1, 4

**Legend: Frequency:** 1. Before direct patient contact 2. After direct patient contact 3. Before and after direct patient contact 4. End of every day. **Disinfection Agent:** 1. 70% Isopropyl alcohol 2. 6% Hydrogen peroxide 3. CL MPS. 4. MP Surface cleaning agent **Who:** A. Optometrist B. Staff.